

TALON® II Litter Series

Talon® II 90C Collapsible Handle Litter: PN:60-0002 NSN:6530-01-504-9051 C€
T3™ Lightweight Assault Litter: PN:60-0022 C€

INSTRUCTIONS FOR USE

● STEP 1:

While standing upright, release buckles from litter.



● STEP 2:

Place litter on the ground and completely extend with the fabric side up.

NOTE:

When litter is in the relaxed position, the hinges closest to you (when holding the handles) will look like Figure 1.



Figure 1



Relaxed Position

EC	REP
----	-----

Manufactured in the
United States of America

REV:030010

PRODUCTS WITH A MISSION®



NORTH AMERICAN RESCUE®


www.NARescue.com • Green, SC 29650 • 888.689.6277

STEP 3:

Keeping the litter straight as possible, grab handles and rotate inward until all hinges rotate & lock (see Figure 2).

NOTE:

This action is best executed simultaneously using 2 individuals on each end of the litter.



STEP 4:

While maintaining hinges in the locked position, apply firm, steady pressure on the spreader bar with foot. Increase pressure with foot until spreader bar locks into place.

NOTE:

To close and store litter, repeat the directions in reverse order.



Collapsible Handle Option:

Stand upright and extend handles to full length.

NOTE:

The button will pop out when the handle is securely extended. Depress the button in order to collapse the handles.



NAR Litter



MAINTENANCE INSPECTION

Although all NAR litters are designed for use in rugged conditions, they still may require regular preventative maintenance.

STEPS FOR CLEANING

- Wipe off all surface dirt with a sponge slightly dampened with plain water
- For more resistant dirt, use a mild solution of water and commercial soap or detergent
- Work up a thick lather using a vigorous back and forth motion
- Rinse thoroughly to remove detergent
- Wipe each belt with a clean cloth and use a mild detergent solution to remove stains
- Hang belts freely to air dry but avoid excessive heat
- Do not store device while still damp
- Storage areas should be clean, dry and free from exposure to fumes or corrosive elements

PRODUCTS WITH A MISSION[®]



North American Rescue, Inc.
10000 Highway 101
Greer, SC 29615

PN:ZZ-0234 REV:022310



NORTH AMERICAN RESCUE[®]

www.NARescue.com • Greer, SC 29650 • 888.639.6277



SAFETY INSPECTION

To maintain service life and high performance, all belts, hinges, handles, and other materials must be inspected frequently. Periodic tests are to be conducted by a trained inspector. Inspect your equipment daily and replace it if any defective conditions are found.

- Examine each belt. Beginning at one end, holding the body side of the belt toward you, grasp the belt with your hands 6 to 8 inches apart.
- Bend the belt into an inverted U shape. This procedure will allow you to locate damaged fibers or cuts more readily. Watch for frayed edges, broken fibers, pulled stitches, cuts, or chemical damage.
- Special attention should be given to the attachment of buckles where attached to webbing. Note any unusual wear, frayed or cut fibers, or distortion of buckles.
- Look for frayed or broken fiber strands that appear as tufts of fiber on the webbing surface.
- Inspect frame, bars, and welds points for cracks, distortions, or stress marks.
- Check pins to ensure they are in place and have fasteners attached properly.
- Inspect handles for cracks, distortions, or stress marks.
- Examine litter bed for general rips, tears, and thread damage. Also look for signs of damage caused by heat, chemicals, or molten metal.

PRODUCTS WITH A MISSION®



NORTH AMERICAN RESCUE®

www.NARescue.com • Greer, SC 29650 • 888.699.6277

PRODUCTS WITH A MISSION

The 3X Standard™ designates absolute "best-of-breed" Combat Casualty Care products. Specifically designed and conceived to sustain the warfighter's life from point of injury to the surgical table. Each item has been manufactured to the most demanding specifications and field-tested in the harshest environmental conditions. Each is a "Product with a Mission", bearing our unconditional guarantee for unsurpassed workmanship and performance.

Our reputation rests on these products. They define who we are and differentiate us from other manufacturers of combat casualty care equipment.

THE 3X STANDARD
 EXACT • EXTREME • EXTRAORDINARY

Visit www.NARescue.com for more info.



North American Rescue takes great pride that this bag is entirely comprised of U.S. materials, manufacturing, and labor and is compliant with the Berry Amendment.

Proudly Made In America


MADE IN AMERICA



STANDARD
E
EXACT • EXTREME
EXTRAORDINARY
NORTH AMERICAN RESCUE, LLC

REV012710

NORTH AMERICAN RESCUE
35 Tedwall Court • Green, SC 29650
888.689.6277 • www.NARescue.com





REV012710

NORTH AMERICAN RESCUE
35 Tedwall Court • Green, SC 29650
888.689.6277 • www.NARescue.com



mPale
A North American Rescue® Exclusive Treatment

This North American Rescue® specialty item is protected against microbial contamination by the revolutionary mPale® antimicrobial with AEGIS™ Shield technology:

Engineered to microscopically protect surfaces. Effective control of bacteria, fungi, algae, mold & mildew. Permanently bonded to the surface. Non-leaching, environmentally friendly. Worldwide proven protection since 1978.



NORTH AMERICAN RESCUE
www.NARescue.com • 888.689.6277
REV02072010